

NEDD8

Reactivity:Human

Tested applications:WB IF IP

Observed MW:Refer to figures

Immunogen:

Recombinant protein

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A11032

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:4738

Isotype:IgG

Swiss Prot:Q15843

Purity:Affinity purification

For research use only.