

## HSPD1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF FC

**Observed MW:**Refer to figures

**Immunogen:**

Synthetic Peptide

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11051

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**3329

**Isotype:**IgG

**Swiss Prot:**P10809

**Purity:**Affinity purification

For research use only.