

NRAS

Reactivity:Human Mouse Rat Zebrafish

Tested applications:WB IP FC

Observed MW:Refer to figures

Immunogen:

Recombinant protein

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A11059

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:4893

Isotype:IgG

Swiss Prot:P01111/P01112/P01116

Purity:Affinity purification

For research use only.