

## EZH2

---

**Reactivity:**Human

**Tested applications:**WB IF

**Recommended Dilution:**WB 1:500 - 1:2000 IF 1:20 - 1:100

**Calculated MW:**85kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human EZH2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

WVS; ENX1; EZH1; KMT6; WVS2; ENX-1; EZH2b; KMT6A

**Catalog #:**A11085

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2146

**Isotype:**IgG

**Swiss Prot:**Q15910

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene.

*To place an order, please [Click HERE](#).*