

## PIK3R2

---

**Reactivity:**Human Rat

**Tested applications:**WB IHC IF IP FC

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11086

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**5296

**Isotype:**IgG

**Swiss Prot:**O00459

**Purity:**Affinity purification

For research use only.