

## ALAS1

---

**Reactivity:**Human

**Tested applications:**WB ICC IF

**Recommended Dilution:**WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**71 kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ALAS1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ALAS; MIG4; ALAS3; ALASH; ALAS-H;

**Catalog #:**A11120

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**211

**Isotype:**IgG

**Swiss Prot:**P13196

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes the mitochondrial enzyme which catalyzes the rate-limiting step in heme (iron-protoporphyrin) biosynthesis. The enzyme encoded by this gene is the housekeeping enzyme; a separate gene encodes a form of the enzyme that is specific for erythroid tissue. The level of the mature encoded protein is regulated by heme: high levels of heme down-regulate the mature enzyme in mitochondria while low heme levels up-regulate. A pseudogene of this gene is located on chromosome 12. Alternative splicing results in multiple transcript variants encoding different isoforms.

*To place an order, please [Click HERE](#).*