

## PTPRB

---

**Reactivity:**Human Mouse

**Tested applications:**IHC IF

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11140

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**5787

**Isotype:**IgG

**Swiss Prot:**P23467

**Purity:**Affinity purification

For research use only.