

RAB4A

Reactivity: Human Mouse Rat

Tested applications: WB ICC IF IP

Recommended Dilution: WB 1:500 - 1:2000 ICC 1:50 - 1:200 IF 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 24 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAB4A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAB4; HRES1; HRES-1; HRES-1/RAB4;

Catalog #: A11143

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 5867

Isotype: IgG

Swiss Prot: P20338

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).