

RAB5A

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:200

Calculated MW:24 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RAB5A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAB5;

To place an order, please [Click HERE](#).

Catalog #:A11146

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:5868

Isotype:IgG

Swiss Prot:P20339

Purity:Affinity purification

For research use only.