

IL18

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IL18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ab

Synonym:

IL18;IGIF;IL-18;IL-1g;IL1F4;MGC12320;

Catalog #: A1115

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3606

Isotype: IgG

Swiss Prot: Q14116

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a proinflammatory cytokine that augments natural killer cell activity in spleen cells, and stimulates interferon gamma production in T-helper type I cells.

Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2011]

To place an order, please [Click HERE](#).