

## RELA

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**65kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human RELA

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

p65; NFKB3;

**Catalog #:**A11204

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5970

**Isotype:**IgG

**Swiss Prot:**Q04206

**Purity:**Affinity purification

For research use only.

**Background:**

NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*