

## VTN

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC ICC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**54/75 kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human VTN

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

VN; V75; VNT;

**Background:**

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11245

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**7448

**Isotype:**IgG

**Swiss Prot:**P04004

**Purity:**Affinity purification

For research use only.