

## ASH2L

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC ICC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**69 kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ASH2L

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ASH2; Bre2; ASH2L1; ASH2L2;

*To place an order, please [Click HERE](#).*

**Catalog #:**A11278

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**9070

**Isotype:**IgG

**Swiss Prot:**Q9UBL3

**Purity:**Affinity purification

For research use only.