

TOMM20

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC IF IP FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW:16 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TOMM20

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAS20; MOM19; TOM20;

To place an order, please [Click HERE](#).

Catalog #:A11308

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:9804

Isotype:IgG

Swiss Prot:Q15388

Purity:Affinity purification

For research use only.