

## DNAJC15

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**16 kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human DNAJC15

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MCJ; HSD18; DNAJD1;

*To place an order, please [Click HERE](#).*

**Catalog #:**A11380

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**29103

**Isotype:**IgG

**Swiss Prot:**Q9Y5T4

**Purity:**Affinity purification

For research use only.