

MKI67

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:359kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human MKI67

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KIA; MIB; MIB-1; PPP1R105;

Background:

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

To place an order, please [Click HERE](#).

Catalog #:A11390

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4288

Isotype:IgG

Swiss Prot:P46013

Purity:Affinity purification

For research use only.