

## MKI67

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**359kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human MKI67

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KIA; MIB; MIB-1; PPP1R105;

**Background:**

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11390

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4288

**Isotype:**IgG

**Swiss Prot:**P46013

**Purity:**Affinity purification

For research use only.