

KRT20

Reactivity:Human Rat

Tested applications:WB IHC ICC IP FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IP 1:20 - 1:50
FC 1:10 - 1:50

Calculated MW:44 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KRT20

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

K20; CD20; CK20; CK-20; KRT21;

Catalog #:A11401

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:54474

Isotype:IgG

Swiss Prot:P35900

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This cytokeratin is a major cellular protein of mature enterocytes and goblet cells and is specifically expressed in the gastric and intestinal mucosa. The type I cytokeratin genes are clustered in a region of chromosome 17q12-q21.

To place an order, please [Click HERE](#).