

## NXPE3

---

**Reactivity:**Human

**Tested applications:**WB IF

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**A11434

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**91775

**Isotype:**IgG

**Swiss Prot:**Q969Y0

**Purity:**Affinity purification

For research use only.