

PINK1

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW:63kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PINK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BRPK; PARK6

Catalog #:A11435

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:65018

Isotype:IgG

Swiss Prot:Q9BXM7

Purity:Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from stress-induced mitochondrial dysfunction. Mutations in this gene cause one form of autosomal recessive early-onset Parkinson disease.

To place an order, please [Click HERE](#).