

RBFOX3

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC IF FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:200
FC 1:20 - 1:50

Calculated MW:34 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RBFOX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol,
pH7.3.

Synonym:

FOX3; NEUN; FOX-3; HRNBP3;

To place an order, please [Click HERE](#).

Catalog #:A11449

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:146713

Isotype:IgG

Swiss Prot:A6NFN3

Purity:Affinity purification

For research use only.