

TP53

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TP53

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A11545

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7157

Isotype: IgG

Swiss Prot: P02340

Purity: Affinity purification

For research use only.