

APOBEC3D

Reactivity:Human Mouse Rat

Tested applications:WB

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human APOBEC3D

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A11598

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:140564

Isotype:IgG

Swiss Prot:Q96AK3

Purity:Affinity purification

For research use only.