

## SLC6A5

---

**Tested applications:**WB

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SLC6A5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11647

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9152

**Isotype:**IgG

**Swiss Prot:**Q9Y345

**Purity:**Affinity purification

For research use only.