

SLC6A2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:69kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SLC6A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NET; NAT1; NET1; SLC6A5;

Background:

This gene encodes a member of the sodium:neurotransmitter symporter family. This member is a multi-pass membrane protein, which is responsible for reuptake of norepinephrine into presynaptic nerve terminals and is a regulator of norepinephrine homeostasis. Mutations in this gene cause orthostatic intolerance, a syndrome characterized by lightheadedness, fatigue, altered mentation and syncope. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.

To place an order, please [Click HERE](#).

Catalog #:A1421

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6530

Isotype:IgG

Swiss Prot:P23975

Purity:Affinity purification

For research use only.