

APOA5

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human APOA5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAP3; APOAV;

Catalog #:A1424

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:116519

Isotype:IgG

Swiss Prot:Q6Q788

Purity:Affinity purification

For research use only.

Background:

Apolipoprotein A5 (ApoA5) is fast gaining attention as a key regulator of serum triglyceride concentrations. An ApoA5 mouse knock-out model produced an approximately four fold increase in serum triglycerides, whereas a knock-in model with human ApoA5 produced 5070% lower concentrations of mouse serum triglycerides. In addition, peroxisome proliferator-activated receptor- agonists, which are used clinically to lower serum triglyceride concentrations, cause increased ApoA5 mRNA expression. Despite these compelling molecular biology data, relatively little is known about ApoA5 protein in human serum. This antibody pair detected recombinant apoA5 protein in sandwich ELISA format and could be potential reagents for the development of clinical diagnostic kits.

To place an order, please [Click HERE](#).