

SELP

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:91kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SELP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM;

Catalog #:A1425

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6403

Isotype:IgG

Swiss Prot:Q6IQ21

Purity:Affinity purification

For research use only.

Background:

SELP is a selenoprotein containing multiple selenocysteine (Sec) residues. The 3' UTR of selenoproteins have a common stem-loop structure, the sec insertion sequence (SECIS), which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. SELP is an extracellular glycoprotein, and it is unusual in that it contains 10 Sec residues per polypeptide. It is a heparin-binding protein that appears to be associated with endothelial cells, and it has been implicated to function as an antioxidant in the extracellular space (1).

To place an order, please [Click HERE](#).