

DPP4

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:110kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DPP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD26; ADABP; ADCP2; DPPIV; TP103;

Background:

CD26 is identical to adenosine deaminase complexing protein-2, and to the T-cell activation antigen CD26. It is an intrinsic membrane glycoprotein and a serine exopeptidase that cleaves X-proline dipeptides from the N-terminus of polypeptides (1).

To place an order, please [Click HERE](#).

Catalog #:A1455

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1803

Isotype:IgG

Swiss Prot:P27487

Purity:Affinity purification

For research use only.