

CD28

Reactivity:Human

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:20 - 1:50

Calculated MW:25kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CD28

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

T-cell-specific surface glycoprotein CD28;TP44;CD_antigen=CD28; Flags;Precursor;

Catalog #:A1468

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:940

Isotype:IgG

Swiss Prot:P10747

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).