

FSHR

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:78kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human FSHR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FSHR;FSHRO;LGR1;MGC141667;MGC141668;ODG1 ;

Catalog #:A1480

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2492

Isotype:IgG

Swiss Prot:P23945

Purity:Affinity purification

For research use only.

Background:

FSHR, also named as LGR1, belongs to the G-protein coupled receptor 1 family and FSH/LSH/TSH subfamily. FSHR is a receptor for follicle-stimulating hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase. Defects in FSHR are a cause of ovarian dysgenesis type 1 (ODG1). Defects in FSHR are a cause of ovarian hyperstimulation syndrome (OHSS). The antibody is special to FSHR.

To place an order, please [Click HERE](#).