

## PICK1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**47kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human PICK1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MGC15204; PICK; PRKCABP

**Catalog #:**A1519

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9463

**Isotype:**IgG

**Swiss Prot:**Q9NRD5

**Purity:**Affinity purification

For research use only.

**Background:**

Protein interacting with C kinase 1 (PICK1) is a PDZ-domain containing protein that is located in the perinuclear region and is phosphorylated in response to PKC

*To place an order, please [Click HERE](#).*