

AIF1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:17kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AIF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IBA1; IRT1; AIF-1; IRT-1;

Catalog #:A1527

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:199

Isotype:IgG

Swiss Prot:P55008

Purity:Affinity purification

For research use only.

Background:

This gene is induced by cytokines and interferon. Its protein product is thought to be involved in negative regulation of growth of vascular smooth muscle cells, which contributes to the anti-inflammatory response to vessel wall trauma. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008]

To place an order, please [Click HERE](#).