

PRKAR2A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:43kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PRKAR2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Synonym:

MGC3606; PKR2; PRKAR2;

Catalog #:A1531

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5576

Isotype:IgG

Swiss Prot:P13861

Purity:Affinity purification

For research use only.

Background:

The effects of cAMP in most tissues and cell types are mainly modulated via protein kinase A, a heterotetrameric protein complex consisting of two regulatory

To place an order, please [Click HERE](#).