

PAH

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:52kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PAH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PH; PKU; PKU1

Catalog #:A1559

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5053

Isotype:IgG

Swiss Prot:P00439

Purity:Affinity purification

For research use only.

Background:

PAH encodes the enzyme phenylalanine hydroxylase that is the rate-limiting step in phenylalanine catabolism. Deficiency of this enzyme activity results in the autosomal recessive disorder phenylketonuria.

To place an order, please [Click HERE](#).