

LGALS1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:1000 - 1:4000 IHC 1:50 - 1:200

Calculated MW:15kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LGALS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GBP; GAL1;

Catalog #:A1580

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3956

Isotype:IgG

Swiss Prot:P09382

Purity:Affinity purification

For research use only.

Background:

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.

To place an order, please [Click HERE](#).