

THRB

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:53kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human THRB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

GRTH; PRTH; THR1; ERBA2; NR1A2; THRB1; THRB2; C-ERBA-2; C-ERBA-BETA;

Catalog #:A1582

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7068

Isotype:IgG

Swiss Prot:P10828

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a nuclear hormone receptor for triiodothyronine. It is one of the several receptors for thyroid hormone, and has been shown to mediate the biological activities of thyroid hormone. Knockout studies in mice suggest that the different receptors, while having certain extent of redundancy, may mediate different functions of thyroid hormone. Mutations in this gene are known to be a cause of generalized thyroid hormone resistance (GTHR), a syndrome characterized by goiter and high levels of circulating thyroid hormone (T3-T4), with normal or slightly elevated thyroid stimulating hormone (TSH). Several alternatively spliced transcript variants encoding the same protein have been observed for this gene. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).