

ADORA2A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ADORA2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADORA2; RDC8; hA2Ar;

Background:

This gene encodes a protein which is one of several receptor subtypes for adenosine. The activity of the encoded protein, a G-protein coupled receptor family member, is mediated by G proteins which activate adenylyl cyclase. The encoded protein is abundant in basal ganglia, vasculature and platelets and it is a major target of caffeine. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).

Catalog #: A1587

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 135

Isotype: IgG

Swiss Prot: P29274

Purity: Affinity purification

For research use only.