

IL23R

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:71kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human IL23R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A1613

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:149233

Isotype:IgG

Swiss Prot:Q5VWK5

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a subunit of the receptor for IL23A/IL23. This protein pairs with the receptor molecule IL12RB1/IL12Rbeta1, and both are required for IL23A signaling. This protein associates constitutively with Janus kinase 2 (JAK2), and also binds to transcription activator STAT3 in a ligand-dependent manner. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).