

SELL

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:42kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic Peptide of human SELL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1;

Background:

This gene encodes a cell surface adhesion molecule that belongs to a family of adhesion/homing receptors. The encoded protein contains a C-type lectin-like domain, a calcium-binding epidermal growth factor-like domain, and two short complement-like repeats. The gene product is required for binding and subsequent rolling of leucocytes on endothelial cells, facilitating their migration into secondary lymphoid organs and inflammation sites. Single-nucleotide polymorphisms in this gene have been associated with various diseases including immunoglobulin A nephropathy.

Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Oct 2009]

To place an order, please [Click HERE](#).

Catalog #:A1622

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6402

Isotype:IgG

Swiss Prot:P14151

Purity:Affinity purification

For research use only.