

CSF1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:60kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CSF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MCSF; CSF-1;

Catalog #:A1627

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1435

Isotype:IgG

Swiss Prot:P09603

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).