

GATA3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GATA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

HDR; HDRS;

Catalog #: A1638

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2625

Isotype: IgG

Swiss Prot: P23771

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.

To place an order, please [Click HERE](#).