

OLR1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:31kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human OLR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LOX1; LOXIN; SLOX1; CLEC8A; SCARE1;

Catalog #:A1639

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4973

Isotype:IgG

Swiss Prot:P78380

Purity:Affinity purification

For research use only.

Background:

This gene encodes a low density lipoprotein receptor that belongs to the C-type lectin superfamily.

This gene is regulated through the cyclic AMP signaling pathway. The encoded protein binds, internalizes and degrades oxidized low-density lipoprotein. This protein may be involved in the regulation of Fas-induced apoptosis. This protein may play a role as a scavenger receptor.

Mutations of this gene have been associated with atherosclerosis, risk of myocardial infarction, and may modify the risk of Alzheimer's disease. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).