

CAMP

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:19kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CAMP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LL37; CAP18; CRAMP; HSD26; CAP-18; FALL39; FALL-39;

Background:

This gene encodes a member of an antimicrobial peptide family, characterized by a highly conserved N-terminal signal peptide containing a cathelin domain and a structurally variable cationic antimicrobial peptide, which is produced by extracellular proteolysis from the C-terminus.

The encoded protein has several functions in addition to antimicrobial activity, including cell chemotaxis, immune mediator induction and inflammatory response regulation. [provided by RefSeq, Aug 2011]

To place an order, please [Click HERE](#).

Catalog #:A1640

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:820

Isotype:IgG

Swiss Prot:P49913

Purity:Affinity purification

For research use only.