

DEFB4A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 7kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DEFB4A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BD-2; SAP1; DEFB2; DEFB4; HBD-2; DEFB-2; DEFB102;

Catalog #: A1643

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1673

Isotype: IgG

Swiss Prot: O15263

Purity: Affinity purification

For research use only.

Background:

Defensins form a family of microbicidal and cytotoxic peptides made by neutrophils. Members of the defensin family are highly similar in protein sequence. This gene encodes defensin, beta 4, an antibiotic peptide which is locally regulated by inflammation.

To place an order, please [Click HERE](#).