

HHEX

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HHEX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HEX; PRH; HMPH; PRHX; HOX11L-PEN;

Catalog #: A1663

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3087

Isotype: IgG

Swiss Prot: Q03014

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the homeobox family of transcription factors, many of which are involved in developmental processes. Expression in specific hematopoietic lineages suggests that this protein may play a role in hematopoietic differentiation.

To place an order, please [Click HERE](#).