

TAC1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:15kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TAC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NK2; NPK; NKNA; TAC2; Hs.2563;

Background:

This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).

Catalog #:A1666

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6863

Isotype:IgG

Swiss Prot:P20366

Purity:Affinity purification

For research use only.