

VTN

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human VTN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VN; V75; VNT;

Background:

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).

Catalog #:A1667

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7448

Isotype:IgG

Swiss Prot:P04004

Purity:Affinity purification

For research use only.