

## CD38

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**34kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant Protein of human CD38

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

T10;

**Catalog #:**A1680

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**952

**Isotype:**IgG

**Swiss Prot:**P28907

**Purity:**Affinity purification

For research use only.

**Background:**

CD38 is a novel multifunctional ectoenzyme widely expressed in cells and tissues especially in leukocytes. CD38 also functions in cell adhesion, signal transduction and calcium signaling.

*To place an order, please [Click HERE](#).*