

GTF2B

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GTF2B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TF2B; TFIIB;

Background:

This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.

To place an order, please [Click HERE](#).

Catalog #:A1708

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2959

Isotype:IgG

Swiss Prot:Q00403

Purity:Affinity purification

For research use only.