

APOA2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:11kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human APOA2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

apoAII; Apo-AII; ApoA-II;

Catalog #:A1711

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:336

Isotype:IgG

Swiss Prot:P02652

Purity:Affinity purification

For research use only.

Background:

This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with apolipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia.

To place an order, please [Click HERE](#).